

KAT14 Protein, Human, Recombinant (His)

General Information

Synonyms:	ADA2A-containing complex subunit 2 (ATAC2);KAT14;CSRP2-binding protein;CSRP2BP;Lysine acetyltransferase 14;CRP2-binding partner (CRP2BP);Cysteine-rich protein 2-binding protein
Protein Construction:	1-782 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q9H8E8
Molecular Weight:	94.9 kDa (predicted)
AA Sequence:	MDSSIHLSLISRHDDEATRSTSEGLEEGEVEGETLLIVESEDQASVDLSHDQSGDSLNSDEGDVSWMEEQL SYFCDKCQKWIPASQLREQLSYLKGDNFFRFTCSDCSADGKEQYERLKLWQQVVMLAMYNSLEGSGRQG YFRWKEDICAFIEKHWTFLGNRKKSTWWSTVAGCLSVGSPMYFRSGAQEFGEPEGWWKLVHNKPPTMKPE GEKLSASTLKIKAAASKPTLDPIITVEGLRKRASRNVPESAMELKEKRSRTQEAKDIRRAQKEAAGFLDRSTSSTP VKFISRGRRPDVILEKGEVIDFSSLSSSDRTPLTSPSPSPSLDFSAPGTPASHSATPSLLSEADLIPDVMPPQALF HDDDEMEGDGVIDPGMEYVPPPAGSVASGPVVGVRKKVIRGPEQIKQVESEEEKPDMDIDSEDTSNTSL QTRAREKRKPQLEKDTKPKEPRYPVSIYEEKLLLKRLEACPGAVAMTPEARRLKRKLIVRQAKRDRGLPLFDL DQVVNAALLLVLDGIYGAKEGGISRLPAGQATYRTTCQDFRILDYQTSLSRKGFRHQTKFLYRLVGSEDMA VDQSIVSPYTSRILKPYIRRDYETKPPKLQLLSQIRSHLHRSDPHWTPEDAPLDYCYVRPNHIPTINSMCQEFF WPGIDLSECLQYPDFSVVLYKKVIAAFGMVDPVKYNEAYISFLVHPEWRRAGIATFMIYHLIQTCMGKDVTL HVSASNPAMLLYQKFGFKTEEVLDYDKYYPLESTECKHAFFLRLRR

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Component of the ATAC complex, a complex with histone acetyltransferase activity on histones H3 and H4. May function as a scaffold for the ATAC complex to promote ATAC complex stability. Has also weak histone acetyltransferase activity toward histone H4. Required for the normal progression through G1 and G2/M phases of the cell cycle.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481