

## Compound T130978

### Chemical Properties


CAS No. :

Formula: C<sub>30</sub>H<sub>38</sub>O<sub>17</sub>

Molecular Weight: 670.62

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description

(2S,3R,4R,5S,6S)-3,5-dihydroxy-2-[[[(2S,4S,5S,6R,10S)-2-(hydroxymethyl)-10-[[[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]-3,9-dioxatricyclo[4.4.0.0<sup>2,4</sup>]dec-7-en-5-yl]oxy]-6-methyloxan-4-yl (2E)-3-(3,4-dihydroxyphenyl)prop-2-enoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4912 mL	7.4558 mL	14.9116 mL
5 mM	0.2982 mL	1.4912 mL	2.9823 mL
10 mM	0.1491 mL	0.7456 mL	1.4912 mL
50 mM	0.0298 mL	0.1491 mL	0.2982 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481