

LILRB3 Protein, Mouse, Recombinant (aa 25-642, His)

General Information

Synonyms:	PIR-B;Cell-surface glycoprotein p91;Leukocyte immunoglobulin-like receptor subfamily B member 3 (LIR-3;Leukocyte immunoglobulin-like receptor 3);Lilrb3;Pirb;Paired immunoglobulin-like receptor B
Protein Construction:	25-642 aa
Species:	Mouse
Expression Host:	HEK293 Cells
Accession:	P97484
Molecular Weight:	70.8 kDa (predicted)
AA Sequence:	SLPKPILRVQPDSVVSRRTKVTLCEETIGANEYRLYKDGKLYKTVTKNKQKPENKAEFSSNVDSL NAGQYRC SYSTQYKSSGYSDDLLELVVTGHYWTPLLAQASPVVTSGGYVTLQCESWHNDHKFILTVEGPQKLSWTQDSQ YNYSTRKYHALFSVGPVTPNQRWICRCYSYDRNRPYVWSPPESEVELLVSGNLQKPTIKAEPG SVITSKRAMTI WCQGNLDAEVYFLHNEKSQKTQSTQTLQEPGNKGGKFFIPSVTLQHAGQYRCYCYGSAGWSQPSDTLELVVT GIYEYEPRLSVLPSPVVTAGGNMTLHCASDFPYDKFILTKEDKKFGNSLDTEHISSSGQYRALFIIGPTTPTHTG AFRCYGYKKNAPQLWSVPSALQQLISGLSKKPSLLTHQHGHILDPGMTLTLQCFSDIN YDRFALHKVGGADIM QHSSQQTDTGFSVANFTLGYVSSSTGGQYRCYGAHNLSSEWSASSEPLDILITGQLPLTPSLSVQPNHTVHSG ETVSLLCWSMDSVDTFILSKEGSAQQPLRLKSKSHDQQSQAEFSMSAVTSHLSGTYRCYGAQDSSFYLLSSAS APVELTVSGPIETSTPPPTMSMPLGGLHMYLK

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in sterile deionized water. The product concentration should not be less than 100 μg/mL. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

May act as receptor for class I MHC antigens. Becomes activated upon coligation of LILRB3 and immune receptors, such as FCGR2B and the B-cell receptor. Down-regulates antigen-induced B-cell activation by recruiting phosphatases to its immunoreceptor tyrosine-based inhibitor motifs (ITIM).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481