


## Erysotramidine

## Chemical Properties

CAS No. :	52358-58-4
Formula:	C <sub>19</sub> H <sub>21</sub> NO <sub>4</sub>
Molecular Weight:	327.37
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(+)-Erysotramidine shows potent dosedependant antifeedant activity at concentrations 3100 ppm.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0546 mL	15.2732 mL	30.5465 mL
5 mM	0.6109 mL	3.0546 mL	6.1093 mL
10 mM	0.3055 mL	1.5273 mL	3.0546 mL
50 mM	0.0611 mL	0.3055 mL	0.6109 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481