

Rp-Uridine-5'-O-(1-thiotriphosphate) sodium

Chemical Properties

CAS No. :


Formula: C₉H₁₁N₂O₁₄P₃S₄Na

Molecular Weight: 588.14

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

Rp-Uridine-5'-O-(1-thiotriphosphate) (Rp-UTP- α -S) serves as a nucleotide agonist for purinergic P2Y₂ and P2Y₄ receptors, effectively inducing inositol phosphate accumulation within 1321N1 astrocytoma cells that express either P2Y₂ or P2Y₄ receptors, with EC₅₀ values of 5.4 and 27 μ M, respectively.

Solubility Information

Solubility

H₂O: Soluble
(< 1 mg/ml refers to the product slightly soluble or insoluble)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7003 mL	8.5014 mL	17.0028 mL
5 mM	0.3401 mL	1.7003 mL	3.4006 mL
10 mM	0.170 mL	0.8501 mL	1.7003 mL
50 mM	0.034 mL	0.170 mL	0.3401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481