


## Pomalidomide 4'-PEG5-amine

## Chemical Properties

CAS No. :	2357117-23-6
Formula:	C <sub>25</sub> H <sub>36</sub> N <sub>4</sub> O <sub>9</sub> .HCl
Molecular Weight:	573.04
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



## Biological Description

Description	Pomalidomide 4'-PEG5-amine, a cereblon ligand designed for PROTAC research and development, features an integrated E3 ligase ligand and a PEG5 linker ending in an amine for easy conjugation to target protein ligands. It is among a series of functionalized tool molecules specifically crafted for PROTAC R&D.
Targets(IC50)	Others,E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7451 mL	8.7254 mL	17.4508 mL
5 mM	0.349 mL	1.7451 mL	3.4902 mL
10 mM	0.1745 mL	0.8725 mL	1.7451 mL
50 mM	0.0349 mL	0.1745 mL	0.349 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481