


AW68

Chemical Properties

CAS No. : 2247236-61-7
 Formula: C₂₂H₂₁ClN₆
 Molecular Weight: 404.89
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	AW68 is a potential small molecule BRD4 inhibitor that prevents and treats BRD4-related diseases and can be used in cancer research.
Targets(IC50)	Epigenetic Reader Domain

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4698 mL	12.349 mL	24.6981 mL
5 mM	0.494 mL	2.4698 mL	4.9396 mL
10 mM	0.247 mL	1.2349 mL	2.4698 mL
50 mM	0.0494 mL	0.247 mL	0.494 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Yujun Zhao, et al. Brd4 inhibitor, preparation and application thereof. WO2018188660A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481