


VX-702

Chemical Properties

CAS No. : 745833-23-2
 Formula: C₁₉H₁₂F₄N₄O₂
 Molecular Weight: 404.32
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	VX-702 is a highly specific p38 α MAPK inhibitor, 14-fold higher selectivity for the p38 α than p38 β . VX-702 is a small molecule investigational oral anti-cytokine therapy for treatment of inflammatory diseases, specifically rheumatoid arthritis (RA).
Targets(IC50)	Autophagy,p38 MAPK
In vitro	Administering VX-702 (0.1 mg/kg) twice daily exhibited effects similar to those of methotrexate (0.1 mg/kg). Furthermore, the efficacy of administering VX-702 (5 mg/kg) twice daily was akin to that of a once-daily dose of prednisolone (10 mg/kg). The MI/AAR ratio in the group treated with VX-702 (50 mg/kg) was significantly lower compared to both the control group and the group treated with VX-702 at 5 mg/kg. VX-702 has a half-life of 16-20 hours, a clearance volume of 3.75 L, and a distribution volume of 73 L/kg. It selectively inhibits the activity of p38 MAPK without affecting ERKs and JNKs.
In vivo	VX-702 does not inhibit platelet aggregation induced by p38 MAPK activators. Treatment with VX-702 (1 μ M) can completely or partially suppress p38 activity induced by various platelet activators, including thrombin, AYPGKF, SFLLRN, U46619, and collagen, with an IC50 ranging from 4-20 nM. Additionally, VX-702 dose-dependently inhibits the production of IL-6 (IC50: 59 ng/ml), IL-1 β (IC50: 122 ng/ml), and TNF α (IC50: 99 ng/ml).
Kinase Assay	HDAC enzyme assay in vitro: The deacetylase enzyme assay is based on a homogeneous fluorescence release assay. Purified recombinant HDAC enzymes are incubated with MGCD0103 diluted in various concentrations for 10 minutes in assay buffer [25 mM HEPES (pH 8.0), 137 mM NaCl, 1 mM MgCl ₂ , 2.7 mM KCl] at room temperature. The substrate Boc-Lys(ϵ -Ac)-AMC is added to the reaction for further incubation at 37 °C. The concentration of the substrate and the incubation time varies for different isotypes of HDAC enzymes. A 20-min trypsin incubation at room temperature allows the release of the fluorophore from the deacetylated substrate. The fluorescent signal is detected by fluorometer at excitation of 360 nm, emission of 470 nm, and cutoff at 435 nm.

Solubility Information

Solubility	Ethanol: < 1 mg/mL (insoluble or slightly soluble), H ₂ O: < 1 mg/mL (insoluble or slightly soluble),
------------	---

A DRUG SCREENING EXPERT

Solubility	DMSO: 250 mg/mL (618.32 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (4.95 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4733 mL	12.3664 mL	24.7329 mL
5 mM	0.4947 mL	2.4733 mL	4.9466 mL
10 mM	0.2473 mL	1.2366 mL	2.4733 mL
50 mM	0.0495 mL	0.2473 mL	0.4947 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Kuliopulos A, et al. Thromb Haemost, 2004, 92(6), 1387-1393.
- Braddock M, IDDB Meeting Report, 2005, March 14-15.
- Gill A, IDDB Meeting Report, 2002, March 06-08.
- Bhattacharya K, et al. Circulation, 2003, 108(17), 882.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481