


PM 02734

Chemical Properties

CAS No. : 915713-02-9
 Formula: C77H125F3N14O18
 Molecular Weight: 1591.922
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	PM 02734 is an antineoplastic.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6282 mL	3.1409 mL	6.2817 mL
5 mM	0.1256 mL	0.6282 mL	1.2563 mL
10 mM	0.0628 mL	0.3141 mL	0.6282 mL
50 mM	0.0126 mL	0.0628 mL	0.1256 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Tomillero A, Moral MA. Gateways to clinical trials. Methods Find Exp Clin Pharmacol. 2008 Oct;30(8):643-72. PubMed PMID: 19088949.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481