


Clofeverine (free base)

Chemical Properties

CAS No. : 54340-63-5
 Formula: C₁₆H₁₆ClNO₃
 Molecular Weight: 305.76
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Clofeverine (free base) is a highly selective relaxant for gastrointestinal smooth muscles with a beta-receptor stimulative nature. On isolated ileum, the inhibitory effect of clofeverine for spontaneous motility was comparable to that of isoproterenol and 100 times higher than that of papaverine.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2705 mL	16.3527 mL	32.7054 mL
5 mM	0.6541 mL	3.2705 mL	6.5411 mL
10 mM	0.3271 mL	1.6353 mL	3.2705 mL
50 mM	0.0654 mL	0.3271 mL	0.6541 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481