

Parstatin(mouse)

Chemical Properties

CAS No. : 1065756-01-5

Formula: C189H326N58O57S3

Molecular Weight: 4419.19

MGPRRLIIVALGLSLCGPLLSRVPMSPESERTDATVNP

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cell-permeable peptide cleaved from protease-activated receptor 1 (PAR1) upon receptor activation. Attenuates endothelial cell migration and proliferation (IC50 ~ 20 μ M), and induces cell cycle arrest. Promotes activation of caspase-3 and exhibits pro-apoptotic activity in vitro. Inhibits angiogenesis and exhibits cardioprotective activity in vivo.
Targets(IC50)	Protease-activated Receptor

Solubility Information

Solubility	H2O: 1 mg/mL (0.23 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.2263 mL	1.1314 mL	2.2629 mL
5 mM	0.0453 mL	0.2263 mL	0.4526 mL
10 mM	0.0226 mL	0.1131 mL	0.2263 mL
50 mM	0.0045 mL	0.0226 mL	0.0453 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Zania et al (2009) Parstatin, the cleaved peptide on proteinase-activated receptor 1 activation, is a potent inhibitor of activation. J.Pharmacol.Exp.Ther. 328 378 PMID:

Strande et al (2009) Parstatin: a cryptic peptide involved in cardioprotection after ischaemic and reperfusion injury. Cardiovasc.Res. 83 325 PMID:

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481