

Dodovisone D

Chemical Properties


CAS No. :

Formula: C44H42O16

Molecular Weight: 826.804

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dodovisone D is a useful organic compound for research related to life sciences and the catalog number is T125677.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2095 mL	6.0474 mL	12.0948 mL
5 mM	0.2419 mL	1.2095 mL	2.419 mL
10 mM	0.1209 mL	0.6047 mL	1.2095 mL
50 mM	0.0242 mL	0.1209 mL	0.2419 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481