


3-(hydroxymethyl)cyclopentanone

Chemical Properties

CAS No. :	113681-11-1
Formula:	C ₆ H ₁₀ O ₂
Molecular Weight:	114.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Hydroxymethylcyclopentanone is a precursor to synthesize methyl epijasmonate. 3-Hydroxymethylcyclopentanone compounds, which are useful as intermediates in the preparation of HIV chemokine CCR-5 receptor antagonists.
Targets(IC50)	CCR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.7612 mL	43.8059 mL	87.6117 mL
5 mM	1.7522 mL	8.7612 mL	17.5223 mL
10 mM	0.8761 mL	4.3806 mL	8.7612 mL
50 mM	0.1752 mL	0.8761 mL	1.7522 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481