

## Virus-Like Particles (VLPs) isotype control

### General Information

Expression Host: HEK293 Cells  
Molecular Weight: 56.1 kDa (predicted)

### QC Testing

Biological Activity: VLPs Isotype Control Protein  
Endotoxin: < 1.0 EU/µg of the protein as determined by the LAL method.  
Formulation: Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481