

## Mouse CD8+ T Cell Isolation Kit (negative selection)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

TargetMol's Mouse CD8+ T Cell Isolation Kit (negative selection) provides superparamagnetic beads for the isolation of CD8+ T cells from single-cell suspensions of mouse splenocytes or other tissues using a negative sorting method. The principle is to use different biotin-labelled monoclonal antibodies to label non-target cells (non-CD8+ T cells), which are then cleared by streptavidin-labelled magnetic beads for the purpose of mouse CD8+ T cell sorting.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481