

Mag Beads Streptavidin (300 nm), Sterile

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

The streptavidin (SA)/biotin interaction, known for its high binding affinity ($K = 10^{15}$), has been widely applied in the biological field. TargetMol Mag Beads Streptavidin covalently connects SA to the beads surface through protein-protein coupling method. This enables efficient binding of biotinylated antibodies, nucleic acids, and proteins. The product combines superparamagnetic advantages with unique features including uniform size and regular morphology, facilitating the rapid capture of target molecules and magnetic separation. This product can be used with automated equipment for high-throughput technology.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481