

Protein A/G Immunoprecipitation Kit

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

TargetMol Protein A/G Immunoprecipitation kit enables Protein A/G to cover superparamagnetic microspheres with high density by utilizing bio-nano technology. Protein A/G Immunomagnetic Beads can be used for the immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. The kit is provided with optimized pre-made buffers, providing optimal reaction conditions for IP experiments and enhancing the stability. This product is suitable for antigen samples derived from cell lysates, cell culture supernatants, serum, ascites, and other sources.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481