

4-methylnonanoic acid

Chemical Properties


CAS No. : 45019-28-1

Formula: C₁₀H₂₀O₂

Molecular Weight: 172.26

Keep away from moisture, Store under nitrogen
 Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-methylnonanoic acid is a RIFM fragrance and one of the sources of mutton odor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8052 mL	29.0259 mL	58.0518 mL
5 mM	1.161 mL	5.8052 mL	11.6104 mL
10 mM	0.5805 mL	2.9026 mL	5.8052 mL
50 mM	0.1161 mL	0.5805 mL	1.161 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Filho, C.B., Furlan, J.M., de Menezes, C.R. et al. Sample Preparation Methods for Fatty Acid Analysis in Different Raw Meat Products by GC-FID. Food Anal. Methods 16, 749–758 (2023).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481