

Butyl octanoate

Chemical Properties


CAS No. : 589-75-3

Formula: C₁₂H₂₄O₂

Molecular Weight: 200.32

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Butyl octanoate is a synthetic flavor ester with a unique aroma.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.992 mL	24.9601 mL	49.9201 mL
5 mM	0.9984 mL	4.992 mL	9.984 mL
10 mM	0.4992 mL	2.496 mL	4.992 mL
50 mM	0.0998 mL	0.4992 mL	0.9984 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Kumari M, et al. Biotechnological potential of psychrophilic microorganisms as the source of cold-active enzymes in food processing applications. 3 Biotech. 2021 Nov;11(11):479.

Yang, S., et al. Evaluation of the volatile profiles in pulp of 85 apple cultivars (*Malus domestica*) by HS-SPME combined with GC-MS. Food Measure 15, 4215-4225 (2021).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481