

Ethenyl tetradecanoate

Chemical Properties

CAS No. : 5809-91-6


Formula: C16H30O2

Molecular Weight: 254.41

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ethenyl tetradecanoate (VinylMyristate) can be used to synthesize other compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9307 mL	19.6533 mL	39.3066 mL
5 mM	0.7861 mL	3.9307 mL	7.8613 mL
10 mM	0.3931 mL	1.9653 mL	3.9307 mL
50 mM	0.0786 mL	0.3931 mL	0.7861 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

de Oliveira Romera C, et al. Enzymatic Synthesis of a Diene Ester Monomer Derived from Renewable Resource. Appl Biochem Biotechnol. 2019 Nov;189(3):745-759.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481