

FE-200041

Chemical Properties


CAS No. : 212710-26-4

Formula: C₃₀H₄₄N₈O₄

Molecular Weight: 580.72

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	FE-200041 is a peripheral efficacious kappa opioid agonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.722 mL	8.610 mL	17.220 mL
5 mM	0.3444 mL	1.722 mL	3.444 mL
10 mM	0.1722 mL	0.861 mL	1.722 mL
50 mM	0.0344 mL	0.1722 mL	0.3444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481