


## Orismilast

## Chemical Properties

CAS No. :	1353546-86-7
Formula:	C <sub>19</sub> H <sub>15</sub> Cl <sub>2</sub> F <sub>2</sub> N <sub>2</sub> O <sub>7</sub> S
Molecular Weight:	510.29
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	LEO-32731 is a bioactive chemical.
Targets(IC50)	Others,PDE

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9597 mL	9.7983 mL	19.5967 mL
5 mM	0.3919 mL	1.9597 mL	3.9193 mL
10 mM	0.196 mL	0.9798 mL	1.9597 mL
50 mM	0.0392 mL	0.196 mL	0.3919 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481