

Halfordamine

Chemical Properties


CAS No. :

Formula: C₁₂H₁₃NO₄

Molecular Weight: 235.239

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Halfordamine is a useful organic compound for research related to life sciences and the catalog number is T131401.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.251 mL	21.2549 mL	42.5098 mL
5 mM	0.8502 mL	4.251 mL	8.502 mL
10 mM	0.4251 mL	2.1255 mL	4.251 mL
50 mM	0.085 mL	0.4251 mL	0.8502 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481