

## Protein G Immunoprecipitation Kit

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

TargetMol Protein G Immunomagnetic Beads enables Protein G to cover superparamagnetic microspheres with high density by utilizing bio-nano technology. Protein G Immunomagnetic Beads can be used for the immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. The Protein G Immunoprecipitation Kit comes with prepared buffers to provide a relatively better condition for immunoprecipitation and to enhance the stability. This product is suitable for antigen samples derived from cell lysates, cell culture supernatants, serum, ascites, and other sources.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481