

Protein L Immunoprecipitation Kit

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

TargetMol Protein L Immunomagnetic Beads enables Protein L to cover superparamagnetic microspheres with high density by utilizing bio-nano technology. Protein L Immunomagnetic Beads can be used for the immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. The Protein L Immunoprecipitation Kit comes with prepared buffers to provide a relatively better condition for immunoprecipitation and to enhance the stability. This product is suitable for antigen samples derived from cell lysates, cell culture supernatants, serum, ascites, and other sources.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481