

HA Tag Immunoprecipitation Kit

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TargetMol's HA Tag Immunomagnetic Beads specifically bind to proteins with the HA tag, which can be used for immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. The HA Tag Immunoprecipitation Kit comes with prepared buffers to provide a relatively better condition for immunoprecipitation and to enhance the stability. This product is suitable for antigen samples from cell lysates, cell culture supernatants, serum, ascites, etc.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481