

MYC Tag Immunoprecipitation Kit

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at 4°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|-------------|--|
| Description | TargetMol's MYC Tag Immunomagnetic Beads specifically bind to proteins with the MYC tag, which can be used for immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. The MYC Tag Immunoprecipitation Kit comes with prepared buffers to provide a relatively better condition for immunoprecipitation and to enhance the stability. This product is suitable for antigen samples from cell lysates, cell culture supernatants, serum, ascites, etc. |
|-------------|--|

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481