

## GST Tag Immunoprecipitation Kit

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at 4°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

TargetMol's GST Tag Immunomagnetic Beads specifically bind to proteins with the GST tag, which can be used for immunoprecipitation (IP) of proteins, protein complexes, protein-nucleic acid complexes, and other antigens. The GST Tag Immunoprecipitation Kit comes with prepared buffers to provide a relatively better condition for immunoprecipitation and to enhance the stability. This product is suitable for antigen samples from cell lysates, cell culture supernatants, serum, ascites, etc.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481