

HLA-A*24:02&B2M&MAGE-A3 (IMPKAGLLI) Monomer Protein, Human, MHC (His & Avi),

General Information

Synonyms: CT1.3;MAGE-3;MZ2-D;MZ2D;HLA-A2402;HLA-A*24:02;Melanoma-associated antigen 3

Protein Construction: Gly25-Thr305(HLA-A*24:02), Ile21-Met119(B2M) and IMPKAGLLI peptide

Species: Human

Expression Host: HEK293 Cells

Accession: AAA59600.1(HLA-A*24:02)&P61769(B2M)&IMPKAGLLI

Molecular Weight: The protein has a predicted MW of 50.20 kDa. Due to glycosylation, the protein migrates to 55-60 kDa based on Tris-Bis PAGE result.

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC

Endotoxin: < 1.0 EU/µg of the protein as determined by the LAL method.

Formulation: Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS (pH 7.4). Typically, 8% trehalose is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Melanoma antigen gene A3 (MAGE-A3) is one of the most immunogenic cancer testis antigens and is common in various types of cancers. MAGE-A3 can be considered as a predictor for poor prognosis and an option for vaccine immunotherapy in patients with PCa.

Reference

Khalvandi A, et al. Nuclear overexpression levels of MAG-E-A3 predict poor prognosis in patients with prostate cancer. APMIS. 2021 Jun;129(6):291-303. doi: 10.1111/apm.13132. Epub 2021 Apr 13. PMID: 33743542.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481