

HLA-E*01:03&B2M&Peptide (VMAPRTLVL) Tetramer Protein, Human, MHC (His & Avi)

General Information

Synonyms:	MHC class I antigen E;HLAE;sHLA-E;QA1;MHC HLA-E alpha-2.1;MHC;MHC HLA-E alpha-1
Protein Construction:	Gly25-Ile305(HLA-E*01:03), Ile21-Met119(B2M) and VMAPRTLVL peptide. Tetramer is assembled by biotinylated monomer and streptavidin.
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P13747(HLA-E*01:03)&P61769(B2M)&VMAPRTLVL
Molecular Weight:	The protein has a predicted MW of 258 kDa. Due to glycosylation, the protein migrates to 260-265 kDa under Non reducing (N) condition based on Tris-Bis PAGE result.

QC Testing

Biological Activity:	Human NKG2C&CD94, His Tag immobilized on CM5 Chip can bind Human HLA-E*01:03&B2M&Peptide (VMAPRTLVL) Tetramer, His Tag with an affinity constant of 1.65 μ M as determined in SPR assay (Biacore T200).
Purity:	> 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC
Endotoxin:	< 1.0 EU/ μ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing PBS (pH 7.4). Typically, 8% trehalose is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months.

Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

HLA-E is a nonclassical member of the major histocompatibility complex class I gene locus. HLA-E protein shares a high level of homology with MHC Ia classical proteins: it has similar tertiary structure, associates with β 2-microglobulin, and is able to present peptides to cytotoxic lymphocytes. The main function of HLA-E under normal conditions is to present peptides derived from the leader sequences of classical HLA class I proteins, thus serving

for monitoring of expression of these molecules performed by cytotoxic lymphocytes.

Reference

Kanevskiy L, et al. Dimorphism of HLA-E and its Disease Association. Int J Mol Sci. 2019 Nov 4;20(21):5496. doi: 10.3390/ijms20215496. PMID: 31690066; PMCID: PMC6862560.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481