


## 3-Chlorophenoxyacetic acid

## Chemical Properties

CAS No. :	588-32-9
Formula:	C <sub>8</sub> H <sub>7</sub> ClO <sub>3</sub>
Molecular Weight:	186.59
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Chlorophenoxyacetic acid, with CAS No. 588-32-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 3-Chlorophenoxyacetic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Solubility Information

Solubility	DMSO: 37.32 mg/mL (200.01 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3593 mL	26.7967 mL	53.5934 mL
5 mM	1.0719 mL	5.3593 mL	10.7187 mL
10 mM	0.5359 mL	2.6797 mL	5.3593 mL
50 mM	0.1072 mL	0.5359 mL	1.0719 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481