


DJT06001

Chemical Properties

CAS No. : 1628182-40-0
 Formula: C₂₁H₂₀ClN₃O₅
 Molecular Weight: 461.92
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	DJT06001 is a selective Factor Xa inhibitor, reducing thrombus formation with low risk of bleeding.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1649 mL	10.8244 mL	21.6488 mL
5 mM	0.433 mL	2.1649 mL	4.3298 mL
10 mM	0.2165 mL	1.0824 mL	2.1649 mL
50 mM	0.0433 mL	0.2165 mL	0.433 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481