


1-Stearoyl-2-linoleoyl-sn-glycero-3-phosphocholine

Chemical Properties

CAS No. :	27098-24-4
Formula:	C44H84NO8P
Molecular Weight:	786.11
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Stearoyl-2-linoleoyl-sn-glycero-3-phosphocholine is an endogenous metabolite.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2721 mL	6.3604 mL	12.7209 mL
5 mM	0.2544 mL	1.2721 mL	2.5442 mL
10 mM	0.1272 mL	0.636 mL	1.2721 mL
50 mM	0.0254 mL	0.1272 mL	0.2544 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481