

(R)-Finerenone-D5

Chemical Properties


CAS No. :

Formula: C₂₁H₁₇D₅N₄O₃

Molecular Weight: 383.46

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

(R)-Finerenone-D5 is a deuterated compound of (R)-Finerenone. (R)-Finerenone has a CAS number of 1050477-30-9.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6078 mL	13.0392 mL	26.0783 mL
5 mM	0.5216 mL	2.6078 mL	5.2157 mL
10 mM	0.2608 mL	1.3039 mL	2.6078 mL
50 mM	0.0522 mL	0.2608 mL	0.5216 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Jean-Louis H. Dasseux, et al. Ketone compounds and compositions for cholesterol management and related uses. WO2002030860A2.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481