


## Deferasirox-D4

## Chemical Properties

CAS No. :	1133425-75-8
Formula:	C <sub>21</sub> H <sub>11</sub> D <sub>4</sub> N <sub>3</sub> O <sub>4</sub>
Molecular Weight:	377.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Deferasirox-D4 is a deuterated compound of Deferasirox (T1457).
Targets(IC50)	Ferroptosis,Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6498 mL	13.2489 mL	26.4978 mL
5 mM	0.530 mL	2.6498 mL	5.2996 mL
10 mM	0.265 mL	1.3249 mL	2.6498 mL
50 mM	0.053 mL	0.265 mL	0.530 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Sentepe I, Eraslan G. Pharmacokinetic of sulfaclozine in broiler chickens[J]. Food & Chemical Toxicology An International Journal Published for the British Industrial Biological Research Association, 2010, 48(1):0-451.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481