


Dexamethasone-D5

Chemical Properties

CAS No. :	358731-91-6
Formula:	C ₂₂ H ₂₄ D ₅ FO ₅
Molecular Weight:	397.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dexamethasone-D5 is a deuterated compound of Dexamethasone (TMSM-0902).
Targets(IC50)	Glucocorticoid Receptor, Mitophagy, Antibacterial, Antibiotic, Autophagy, Complement System, SARS-CoV

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5158 mL	12.5789 mL	25.1579 mL
5 mM	0.5032 mL	2.5158 mL	5.0316 mL
10 mM	0.2516 mL	1.2579 mL	2.5158 mL
50 mM	0.0503 mL	0.2516 mL	0.5032 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sorkin, E.M., et al.: Drugs, 29, 208 (1985).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481