

Hydroxyzine dihydrochloride-d8

Chemical Properties


CAS No. :

Formula: C₂₁H₂₁D₈Cl₃N₂O₂

Molecular Weight: 455.88

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hydroxyzine-D8 dihydrochloride is a deuterated compound of Hydroxyzine dihydrochloride (T2200).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1936 mL	10.9678 mL	21.9356 mL
5 mM	0.4387 mL	2.1936 mL	4.3871 mL
10 mM	0.2194 mL	1.0968 mL	2.1936 mL
50 mM	0.0439 mL	0.2194 mL	0.4387 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481