

LBQ657-D4

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	LBQ657-D4 is a deuterated compound of LBQ657. LBQ657 (TQ0106) has a CAS number of 149709-44-4. Sacubitrilat (LBQ657) is an effective inhibitor of active neprilysin (NEP).
-------------	--

## Reference

- Schiering N, et al. Structure of neprilysin in complex with the active metabolite of sacubitril. Sci Rep. 2016 Jun 15; 6:27909.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481