


Nandrolone-D3

Chemical Properties

CAS No. :	120813-22-1
Formula:	C ₁₈ H ₂₃ D ₃ O ₂
Molecular Weight:	277.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Nandrolone-D3 is a deuterated compound of Nandrolone.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6046 mL	18.0232 mL	36.0464 mL
5 mM	0.7209 mL	3.6046 mL	7.2093 mL
10 mM	0.3605 mL	1.8023 mL	3.6046 mL
50 mM	0.0721 mL	0.3605 mL	0.7209 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481