

R-(-)-Tretoquinol-D9 hydrochloride

Chemical Properties


CAS No. :

Formula: C₁₉H₁₅D₉ClNO₅

Molecular Weight: 390.91

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	R-(-)-Tretoquinol-D9 hydrochloride is a deuterated compound of R- (-)-Tretoquinol Hydrochloride.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5581 mL	12.7907 mL	25.5813 mL
5 mM	0.5116 mL	2.5581 mL	5.1163 mL
10 mM	0.2558 mL	1.2791 mL	2.5581 mL
50 mM	0.0512 mL	0.2558 mL	0.5116 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481