

## Sacubitril-D4 calcium salt

## Chemical Properties


CAS No. :

Formula: C<sub>48</sub>H<sub>48</sub>D<sub>8</sub>CaN<sub>2</sub>O<sub>10</sub>

Molecular Weight: 869.09

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

Sacubitril-D4 calcium salt is a deuterated compound of Sacubitril calcium salt (T4200).

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1506 mL	5.7531 mL	11.5063 mL
5 mM	0.2301 mL	1.1506 mL	2.3013 mL
10 mM	0.1151 mL	0.5753 mL	1.1506 mL
50 mM	0.023 mL	0.1151 mL	0.2301 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481