

## Tecovirimat-D4

## Chemical Properties


CAS No. :

Formula: C<sub>19</sub>H<sub>11</sub>D<sub>4</sub>F<sub>3</sub>N<sub>2</sub>O<sub>3</sub>

Molecular Weight: 380.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Tecovirimat-D4 is a deuterated compound of Tecovirimat. Tecovirimat (T17026) has a CAS number of 816458-31-8.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6292 mL	13.1458 mL	26.2916 mL
5 mM	0.5258 mL	2.6292 mL	5.2583 mL
10 mM	0.2629 mL	1.3146 mL	2.6292 mL
50 mM	0.0526 mL	0.2629 mL	0.5258 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481