


Torseמידe-D7

Chemical Properties

CAS No. : 1189375-06-1
 Formula: C₁₆H₁₃D₇N₄O₃S
 Molecular Weight: 355.46
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Torseמידe-D7 is a deuterated compound of Torsemide (T1410).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8133 mL	14.0663 mL	28.1326 mL
5 mM	0.5627 mL	2.8133 mL	5.6265 mL
10 mM	0.2813 mL	1.4066 mL	2.8133 mL
50 mM	0.0563 mL	0.2813 mL	0.5627 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481