


## Nebidrazine

## Chemical Properties

CAS No. :	55248-23-2
Formula:	C <sub>9</sub> H <sub>8</sub> Cl <sub>2</sub> N <sub>6</sub>
Molecular Weight:	271.11
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Nebidrazine is an alpha-adrenoceptor agonist.
Targets(IC50)	Others,Histamine Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6885 mL	18.4427 mL	36.8854 mL
5 mM	0.7377 mL	3.6885 mL	7.3771 mL
10 mM	0.3689 mL	1.8443 mL	3.6885 mL
50 mM	0.0738 mL	0.3689 mL	0.7377 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481