


## 1,2-Dilinoleoyl-sn-glycero-3-PC

## Chemical Properties

CAS No. :	6542-05-8
Formula:	C44H80NO8P
Molecular Weight:	782.097
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1,2-Dilinoleoyl-sn-glycero-3-PC (Dilinoleoyllecithin), a phospholipid, finds application in the production of artificial membranes [1].
Targets(IC50)	Others,Liposome

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2786 mL	6.393 mL	12.7861 mL
5 mM	0.2557 mL	1.2786 mL	2.5572 mL
10 mM	0.1279 mL	0.6393 mL	1.2786 mL
50 mM	0.0256 mL	0.1279 mL	0.2557 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481