


1,2,3-Tri- α -Linolenoyl Glycerol

Chemical Properties

CAS No. :	14465-68-0
Formula:	C57H92O6
Molecular Weight:	873.3
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,2,3-Tri- α -linolenoyl glycerol, a polyunsaturated triacylglycerol featuring α -linolenoyl (18:3) side chains at the C-1, C-2, and C-3 positions, serves a key role in lipid research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1451 mL	5.7254 mL	11.4508 mL
5 mM	0.229 mL	1.1451 mL	2.2902 mL
10 mM	0.1145 mL	0.5725 mL	1.1451 mL
50 mM	0.0229 mL	0.1145 mL	0.229 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481