


cis-3,4',5-Trimethoxy-3'-hydroxystilbene

Chemical Properties

CAS No. :	586410-08-4
Formula:	C17H18O4
Molecular Weight:	286.327
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cis-3,4',5-Trimethoxy-3'-hydroxystilbene, a stilbene derivative, induces apoptosis through the mitochondrial release of cytochrome c and suppresses tubulin polymerization. It is also noted for its application in leukemic research [1].
Targets(IC50)	Apoptosis,Others,Microtubule Associated

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4925 mL	17.4624 mL	34.9247 mL
5 mM	0.6985 mL	3.4925 mL	6.9849 mL
10 mM	0.3492 mL	1.7462 mL	3.4925 mL
50 mM	0.0698 mL	0.3492 mL	0.6985 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481