


## 13,14-dihydro-15(R)-Prostaglandin E1

### Chemical Properties

CAS No. : 201848-10-4  
 Formula: C<sub>20</sub>H<sub>36</sub>O<sub>5</sub>  
 Molecular Weight: 356.503  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



### Biological Description

Description	13,14-Dihydro-15(R)-Prostaglandin E1 (13,14-dihydro-15(R)-PGE1) is an analog of 13,14-dihydro-PGE1 characterized by its R-configured hydroxyl group at the C-15 position [1].
Targets(IC50)	Others,Endogenous Metabolite

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.805 mL	14.0252 mL	28.0505 mL
5 mM	0.561 mL	2.805 mL	5.6101 mL
10 mM	0.2805 mL	1.4025 mL	2.805 mL
50 mM	0.0561 mL	0.2805 mL	0.561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481