


13(S)-HOTrE

Chemical Properties

CAS No. :	87984-82-5
Formula:	C ₁₈ H ₃₀ O ₃
Molecular Weight:	294.4
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	13(S)-HOTrE, the 15-lipoxygenase (15-LO) metabolite of linolenic acid, is present in cell membranes and as a cholesteryl ester in atherosclerotic lesions, as well as in the biomembranes of soybeans subjected to 15-LO.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3967 mL	16.9837 mL	33.9674 mL
5 mM	0.6793 mL	3.3967 mL	6.7935 mL
10 mM	0.3397 mL	1.6984 mL	3.3967 mL
50 mM	0.0679 mL	0.3397 mL	0.6793 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481