

MLS001006105

Chemical Properties


CAS No. : 455310-36-8

Formula: C₁₅H₁₅N₅O₃S₂

Molecular Weight: 377.441

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	MLS001006105 is a chemical compound employed in biological assays targeting the protein lamin isoform A-delta10 [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6494 mL	13.2471 mL	26.4943 mL
5 mM	0.5299 mL	2.6494 mL	5.2989 mL
10 mM	0.2649 mL	1.3247 mL	2.6494 mL
50 mM	0.053 mL	0.2649 mL	0.5299 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481